



For the Life Scientist

## AppTaq POLYMERASE

**CONCENTRATION:** 5U/ $\mu$ l

**Store at -20°C.** (The kit will retain full activity for 12 months at -20°C. Can be stored at 4°C for 1 month and go through 30 freeze/thaw cycles with no loss of activity. Avoid prolonged exposure to light).

### DESCRIPTION

AppTaq Polymerase is a highly thermostable enzyme for routine PCR amplification of DNA fragments up to 6 kb. It comes with a 5x reaction buffer which has pre-added enhancers, stabilisers, MgCl<sub>2</sub> and dNTPs to maximise PCR yields, and works in fast or standard thermal cycling conditions. The enzyme generates 3' adenine overhangs on the PCR products which can then be cloned into TA vectors.

### ORDERING INFORMATION

Component	ARP001	ARP003
AppTaq Polymerase (5U/ $\mu$ l)	500 units (1x 0.1ml)	2000 units (4x 0.1ml)
5x Appleton reaction buffer	4x 1ml	16 x 1ml

### PROTOCOL

Prepare a PCR master mix by mixing molecular biology grade water, 5x Appleton reaction buffer, forward and reverse primers and AppTaq Polymerase. Prepare sufficient master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes / wells and then add template DNA.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Place a thin-walled PCR tube/plate on ice and add the following components for each 50 $\mu$ L reaction:

Reagent	Final Concentration	50 $\mu$ L reaction
5x Appleton reaction buffer	1X	10 $\mu$ L
Forward primer (10 $\mu$ M)	400nM	2.0 $\mu$ L
Reverse primer (10 $\mu$ M)	400nM	2.0 $\mu$ L
Template DNA	100 - 500ng*	variable
AppTaq Polymerase (5U/ $\mu$ L)	1.25U - 5U	0.25 $\mu$ L - 1 $\mu$ L
Molecular Biology Grade water, (BMW001)		Up to 50 $\mu$ l final volume

3. Gently vortex the samples and spin down.
4. If using a thermal cycler that does not use a heated lid, overlay the reaction mixture with 25 $\mu$ L of mineral oil.



For the Life Scientist

5. Perform PCR using recommended thermal cycling conditions:

Step	Temperature / °C	Time	Number of cycles
Initial denaturation	95	1-3 min	1
Denaturation	95	15 s	25-40
Annealing	55-65	15 s	
Extension	72	1-90 s (15 s per kb)	

## CONSIDERATIONS

### Template DNA\*

For optimal results, use between 5ng and 500ng per reaction for eukaryotic DNA, and for cDNA use below 100ng in the 50 µL reaction volume. Higher amount of template increases the risk of non-specific PCR products. Trace amounts of certain agents used for DNA purification, such as phenol, EDTA and proteinase K, can inhibit DNA polymerases. Ethanol precipitation and repeated washes of the DNA pellet with 70% ethanol normally removes trace contaminants from DNA samples.

### Reaction Buffer

The 5x Appleton reaction buffer contains optimal concentrations of MgCl<sub>2</sub> (15mM) and dNTPs (5mM), as well as enhancers and stabilisers. This avoids having to vary these parameters to obtain maximum PCR yields. We do not recommend adding further enhancers or magnesium.

### Primers

The recommended concentration range of the PCR primers is 0.1-1 µM. Excessive primer concentrations increase the probability of mis-priming and non-specific PCR products.

### Denaturation

Complete initial denaturation of the template DNA is essential for efficient utilization of the template during the first amplification cycle. If the GC content of the

template is 50% or less, an initial 1-3 min denaturation at 95°C is sufficient. For GC-rich templates, we recommend using Mega-App Polymerase (ARP031) which has been developed specifically for the amplification of extremely difficult templates.

### Annealing

The optimal annealing temperature is 5°C lower than the melting temperature (T<sub>m</sub>) of the primers. Incubation for 0.25-2 min is usually sufficient. However, if non-specific PCR products are obtained in addition to the expected product, the annealing temperature should be optimized by increasing it stepwise by 1-2°C.

### Extension

The optimal extension temperature for AppTaq Polymerase is 70-75°C. The recommended extension step is 30 s at 72°C for PCR products up to 2 kb. For larger products, the extension time should be prolonged by 15 s / kb.

## TROUBLE SHOOTING / TECHNICAL SUPPORT

For troubleshooting please visit [www.appletonwoods.co.uk/PCRtroubleshooting.pdf](http://www.appletonwoods.co.uk/PCRtroubleshooting.pdf) for a trouble shooting guide on PCR. If this does not resolve your issues, please email [technical@appletonwoods.co.uk](mailto:technical@appletonwoods.co.uk) with details of your: amplicon size, reaction setup, cycling conditions, gel images.

**Notes:** AppTaq Polymerase has an error rate of 1 error per 2.0 x 10<sup>5</sup> nucleotides incorporated. AppTaq Polymerase is for research use only.

## ASSOCIATED PRODUCTS

Product	Pack Size	Product Code
Molecular Biology Grade Agarose	500g	AG001
AxyPrep Mag PCR clean up Kit	50mL, 1110 preps	AX402
Molecular biology grade water	500mL	BMW002

More pack sizes available at [www.appletonwoods.co.uk](http://www.appletonwoods.co.uk)